


# Quality Assurance Documentation

<b>ProjectID:</b>	2024AATESD	<b>ConstructID:</b>	24AB3XED
<b>Gene Name:</b>	AzamiRed10	<b>Gene Size:</b>	708 bp
<b>Designation:</b>	E.coli K12 DH10B™ T1R	<b>Vector Backbone:</b>	pMK-RQ (KanR)
<b>Issue Date:</b>	21 June 2024	<b>Quantity:</b>	~5 µg Plasmid DNA

**Product Description:** The synthetic gene AzamiRed10 was assembled from synthetic oligonucleotides and/or PCR products. The fragment was inserted into pMK-RQ (KanR). The plasmid DNA was purified from transformed bacteria and concentration determined by UV spectroscopy. The final construct was verified by sequencing. The sequence identity within the insertion sites was 100%. 5 µg of the plasmid preparation were vacuum dried for shipping.

**Product Handling:** The delivered DNA amounts are indicated on the individual tube labels. Centrifuge tubes prior to opening. Do not store vacuum dried DNA for a prolonged time. Add an appropriate amount of distilled water or 10 mM Tris-HCl (pH 8.5) and incubate for 1 hour at room temperature (optionally followed by an overnight incubation at 4°C). Resuspend DNA by gently pipetting up and down a couple of times. If not to be used immediately, resuspended DNA should be stored at -20°C or -80°C. Storing as aliquots helps to reduce unfavorable freeze-thaw cycles. We recommend sequence verification after each subcloning respectively transformation step.

## Plasmid Map:



## Data Handling:



If you have subscribed to our convenient service to deliver all Quality Assurance Documentation (QAD) for your GeneArt orders electronically through Thermo Fisher Cloud, please follow this link to access your information:

<https://apps.thermofisher.com/apps/geneart-qa>. Alternatively, you can scan the QR code.

If you have not subscribed to this service, the QAD can be found on the CD included in your shipment. The QAD includes a .gb file containing the full vector information, a .fas file for your insert sequence only and if applicable, your ordered raw sequencing data.